


Microfocusing – Recent development for a better use of a microbeam

David FLOT (ESRF)

Microbeams improve signal to noise by minimising solvent diffraction.


(Beam size; 100X100 μm vs 10 μm \varnothing)



- $I(\text{HKL}) \sim I_{\text{beam}} (V_{\text{crystal}} / V_{\text{cell}}) \times (\text{some factors}) \times |F(\text{HKL})|^2$
- If the crystal size is small, if the unit cell is large, we need to increase dramatically the number of incoming photons.
- ➔ we need a small and intense X-ray beam ...

How to achieve a small and intense beam?

How to achieve a small and intense beam?



$$b = S \times (d / D)$$

D = source to focusing element distance, d = focusing element to focus point (sample) distance

S = source size, b = beam size

- **reduce S** (property of the electron beam ...),
- **increase the “source to focusing element distance”** (more sensitive to the angular variation of the incoming beam),
- **reduce the “focusing element to focus point (sample) distance”** (increase of the divergence)


The quality of the focusing device has a contribution to the beam shape.

If we used a mirror, it is the mirror surface (“roughness”) and the mirror shape (“slope error”) which have a contribution to the beam shape.

2 approaches (at the ESRF):

- a moderately focussed beam ($\sim 50 \mu\text{m}$ Ø) + beam defining aperture
- 2 smaller flat mirrors mechanically bended.

1st approach : ID29 (MD2).




-Full beam

- 50 μm \varnothing beam defining aperture


-15 μm \varnothing beam defining aperture

2nd approach : ID23-2 (KB system).





David Flot, ESRF MX School 2010, 9th of
February 2010

ID23-2



Horizontal beamsize: $\sim 7 \mu\text{m}$ (FWHM)
Vertical Beamsize: $< 5 \mu\text{m}$ (FWHM)







Measurement : absorption scan with a W wire (200 μm diameter), 1 μm step

Crystal centring

- Very important (but not always necessary...),
- aim: the crystal should rotate around itself (and not around the beam...),
- Problem: we need to see the crystal (importance of the sample preparation),
- Tactic: align the loop “edge on”, move back by 90° and start the “3 click procedure”:
 - First click: face on (good view),
 - Second click: middle of the loop, on the same vertical line,
 - third click: face on again, on the same vertical line.



Easy case




David Flot, ESRF MX School 2010, 9th of
February 2010



David Flot, ESRF MX School 2010, 9th of
February 2010




David Flot, ESRF MX School 2010, 9th of
February 2010





David Flot, ESRF MX School 2010, 9th of
February 2010

Line scan




Mesh scan for crystal search



5 x 5 , 80 by 80 μm^2 grid

David Flot, ESRF MX School 2010, 9th of
February 2010

1 exemple : the use of the spatial resolution provided by a small beam (Carlo Petosa, IBS-Grenoble)



Sample: domain BD1 (res.17-136) + 20-mer peptide

Crystal ($\sim 30 \times 30 \mu\text{m}^2$) cryoprotected in 30% glycerol


Data Collection: 2/3 bunch 160 mA

Crystals grown from robot screen.


David Flot, ESRF MX School 2010, 9th of February 2010

Images obtained by shooting centre of xtal

Phi=60-61





Phi=150-151





too many spots!

Exposure: 4 sec/image
10% transmission

David Flot, ESRF MX School 2010, 9th of
February 2010



Images obtained by shooting left tip of crystal (volume 1)




Nice-looking spots visible to 2.3

One major (and one minor) lattice visible


Strategy:

Shoot 6 different corners/edges of crystal
Collect 10-30 images per volume
Find combination of images giving best data



Radiation Damage

Image 1
Phi=93-94



spots visible to 2.3






Image 10
Phi=102-103



spots visible to 2.6




Best 90 Images

Completeness: 96.7%

Rfactor : 14.3%

I/Sig (2.4-2.5Å): 2.8




				XSCALE statistics							
RESOLUTION	NUMBER OF REFLECTIONS			COMPLETENESS	R-FACTOR	R-FACTOR	COMPARED	I/SIGMA	R-meas	Rmrgd-F	Mult
LIMIT	OBSERVED	UNIQUE	POSSIBLE	OF DATA	observed	expected					
20.00	46	20	29	69.0%	11.8%	8.9%	46	5.14	14.2%	25.3	2.300
12.00	231	83	93	89.2%	6.0%	7.2%	231	13.16	7.2%	4.5	2.783
8.00	697	245	259	94.6%	7.1%	7.4%	697	14.39	8.4%	5.2	2.845
6.00	1403	456	479	95.2%	9.1%	8.9%	1403	11.89	10.7%	7.1	3.077
5.00	1818	560	579	96.7%	10.2%	9.5%	1818	11.70	12.0%	7.7	3.246
4.50	1604	504	515	97.9%	9.6%	9.2%	1604	12.44	11.3%	7.2	3.183
4.00	2586	770	787	97.8%	9.1%	9.2%	2586	12.53	10.7%	6.9	3.358
3.50	4321	1279	1290	99.1%	10.5%	10.5%	4321	10.50	12.4%	9.4	3.378
3.20	4138	1196	1205	99.3%	13.1%	13.7%	4138	8.18	15.5%	12.9	3.460
3.00	3785	1089	1096	99.4%	17.7%	17.9%	3785	6.49	21.0%	18.6	3.476
2.80	4831	1379	1385	99.6%	21.5%	22.2%	4831	5.39	25.5%	23.9	3.503
2.60	6430	1851	1860	99.5%	29.3%	29.8%	6430	4.19	34.8%	34.1	3.474
2.50	4065	1163	1175	99.0%	37.7%	37.2%	4065	3.48	44.9%	48.6	3.495
2.40	3477	1266	1363	92.9%	41.2%	40.5%	3477	2.82	49.7%	51.4	2.746
2.37	439	281	444	63.3%	41.2%	41.2%	439	2.31	51.5%	51.0	1.562
total	39871	12142	12559	96.7%	14.3%	14.3%	39871	7.11	16.9%	18.4	3.284

- Crystal quality not the same along the crystal,
- A beam smaller than the crystal allows to choose the best place and / or to collect several small datasets,
- “hot beam”: the beam damage is important (no more that 30s total exposure time per point) (*but difficult to predict if the cross section between the crystal and the beam change*).

Small beam / large crystal

- **One crystal, multiple datasets** : available through mxCuBE on all the stations.
- **One crystal, one dataset, data collection between two centred positions**: available on request (will be implemented on a GUI soon).

F AE crystal




Statistics from XDS (Correct step)

SUBSET OF INTENSITY DATA WITH SIGNAL/NOISE >= -3.0 AS FUNCTION OF RESOLUTION													
RESOLUTION	NUMBER OF REFLECTIONS			COMPLETENESS	R-FACTOR	R-FACTOR COMPARED	I/SIGMA	R-meas	Rmrgd-F	Anomal	SigAno	Nano	
LIMIT	OBSERVED	UNIQUE	POSSIBLE	OF DATA	observed	expected				Corr			
5.07	36929	6490	6526	99.4%	2.9%	3.5%	36928	43.82	3.2%	1.6%	90%	3.133	2853
3.61	67454	11654	11654	100.0%	3.1%	3.6%	67454	42.72	3.4%	1.7%	78%	2.172	5425
2.95	87304	15019	15019	100.0%	3.7%	4.1%	87304	36.28	4.1%	2.2%	73%	2.067	7113
2.56	103798	17801	17802	100.0%	4.9%	5.1%	103798	28.66	5.3%	3.1%	69%	1.926	8502
2.29	117267	20173	20173	100.0%	6.1%	6.3%	117267	24.02	6.7%	4.1%	59%	1.637	9687
2.09	128809	22249	22249	100.0%	7.7%	7.7%	128809	19.93	8.5%	5.4%	49%	1.418	10722
1.94	139178	24169	24169	100.0%	10.8%	10.7%	139178	15.15	11.9%	7.8%	38%	1.245	11689
1.81	149092	26055	26055	100.0%	16.0%	16.0%	149092	10.73	17.6%	11.8%	26%	1.075	12622
1.71	106265	25577	27669	92.4%	23.2%	23.4%	103797	6.12	26.3%	22.6%	14%	0.927	10900
total	936096	169187	171316	98.8%	5.6%	5.9%	933627	21.09	6.1%	5.6%	53%	1.523	79513

First data set: normal data collection , 270 degrees, 1 s, 1 degree per image

FAE crystal



First data set: normal data collection , 270 degrees, **1 s** 1 degree per image

Second dataset: “helical” data collection, **1.5 s** exposure time per image

Statistics from XDS (Correct step)


SUBSET OF INTENSITY DATA WITH SIGNAL/NOISE >= -3.0 AS FUNCTION OF RESOLUTION													
RESOLUTION	NUMBER OF REFLECTIONS			COMPLETENESS	R-FACTOR	R-FACTOR COMPARED	I/SIGMA	R-meas	Rmrgd-F	Anomal	SigAno	Nano	
LIMIT	OBSERVED	UNIQUE	POSSIBLE	OF DATA	observed	expected				Corr			
5.07	36929	6490	6526	99.4%	2.9%	3.5%	36928	43.82	3.2%	1.6%	90%	3.133	2853
3.61	67454	11654	11654	100.0%	3.1%	3.6%	67454	42.72	3.4%	1.7%	78%	2.172	5425
2.95	87304	15019	15019	100.0%	3.7%	4.1%	87304	36.28	4.1%	2.2%	73%	2.067	7113
2.56	103798	17801	17802	100.0%	4.9%	5.1%	103798	28.66	5.3%	3.1%	69%	1.926	8502
2.29	117267	20173	20173	100.0%	6.1%	6.3%	117267	24.02	6.7%	4.1%	59%	1.637	9687
2.09	128809	22249	22249	100.0%	7.7%	7.7%	128809	19.93	8.5%	5.4%	49%	1.418	10722
1.94	139178	24169	24169	100.0%	10.8%	10.7%	139178	15.15	11.9%	7.8%	38%	1.245	11689
1.81	149092	26055	26055	100.0%	16.0%	16.0%	149092	10.73	17.6%	11.8%	26%	1.075	12622
1.71	106265	25577	27669	92.4%	23.2%	23.4%	103797	6.12	26.3%	22.6%	14%	0.927	10900
total	936096	169187	171316	98.8%	5.6%	5.9%	933627	21.09	6.1%	5.6%	53%	1.523	79513

SUBSET OF INTENSITY DATA WITH SIGNAL/NOISE >= -3.0 AS FUNCTION OF RESOLUTION													
RESOLUTION	NUMBER OF REFLECTIONS			COMPLETENESS	R-FACTOR	R-FACTOR COMPARED	I/SIGMA	R-meas	Rmrgd-F	Anomal	SigAno	Nano	
LIMIT	OBSERVED	UNIQUE	POSSIBLE	OF DATA	observed	expected				Corr			
5.09	48730	6401	6439	99.4%	2.0%	2.5%	48727	74.62	2.2%	0.9%	97%	5.649	2814
3.61	88981	11607	11607	100.0%	2.2%	2.6%	88981	72.69	2.4%	1.0%	93%	3.836	5403
2.95	115505	15025	15025	100.0%	2.6%	2.8%	115505	63.40	2.8%	1.2%	91%	3.700	7114
2.56	137168	17774	17774	100.0%	3.1%	3.3%	137168	52.18	3.4%	1.6%	89%	3.498	8489
2.29	155443	20164	20164	100.0%	3.7%	3.8%	155443	45.10	4.0%	2.0%	84%	2.963	9680
2.09	171089	22225	22225	100.0%	4.5%	4.5%	171089	38.57	4.8%	2.5%	79%	2.534	10710
1.93	185579	24196	24196	100.0%	5.8%	5.9%	185579	30.71	6.3%	3.4%	72%	2.139	11703
1.81	198774	26041	26041	100.0%	7.9%	8.0%	198774	23.18	8.5%	4.7%	63%	1.762	12616
1.71	141547	25989	27679	93.9%	10.6%	11.0%	139734	13.94	11.7%	8.6%	47%	1.334	11683
total	1242816	169422	171150	99.0%	3.5%	3.7%	1241000	39.41	3.7%	2.6%	79%	2.634	10212

First data set: normal data collection , 270 degrees, 1 s 1 degree per image

Seconda dataset: helical data collection, 1.5 s exposure time per image

B factor from mosflm/scala *



But...

- Crystal centring...
- Crystal homogeneity...
- Cell parameters...

Good practices with a small beam

1. Use a loop which matches the crystal size,
2. Minimise the amount of cryoprotectant, use one which is not too viscous,
3. Come to the station with a diffraction plan,
4. Do a quick centring, move the crystal to the beam, collect one diffraction pattern,
5. Centre the crystal properly if the diffraction quality is good enough,
6. Be care full with **radiation damage** when collecting (no more that 30s total exposure time per point) (*but difficult to predict if the cross section between the crystal and the beam change*),
7. Check the data quality when collecting,
8. When collecting small datasets, use “overlap” between the datasets.

Summary

- Possible to collect data on “**desperate**” **subjects** (balance between more time in the wet lab or more time in the beam line / computer room),
- Small crystals can be **less** affected by **twinning**,
- **Special shapes like needle** like crystals can profit from a small beam,
- **Micro-volumes** of larger crystals can be **better ordered** than the **macrocrystal** (due to twinning, high mosaicity, etc...).

**Thank You
for your attention!**